



Pair of Permeability End Caps NX supplied with Distance Block

Code: [70-1750](#)

Product Group: [Permeability of Rock](#)

Permeability end caps used with ELE-Hoek Cells and constant pressure systems are a cost-effective solution suited to investigating the permeability of rock at high confining pressures in the laboratory. To collect and measure the water which permeates through the rock specimen, a suitable burette such as 25-4540 is recommended. Each permeability end cap incorporates a tubing connector which accepts standard 6 mm tubing, 26-1926, used to connect the cell to the pressure system and burette.

Accessories



Nylon Tubing 6 mm OD x 4 mm ID 3500 Kpa

Code: [26-1926](#)

Alternatives



ELE-Hoek Cell ADR Touch 2000 Compression Machine

Code: [70-2630/01](#)



ELE-Hoek Cell NX 54.74 mm Diameter

Code: [70-1710](#)



Hand Operated Pressure System with Pressure Gauge & Flexible Hose

Code: 70-5000



Pressure Test 3500 Oil/Water Constant Pressure System with Digital Pressure Gauge

Code: 70-5130/01